


BI+ BL6

Wireless Electric Derma Pen



www.portable-ultrasound-scanner.com

Quality, Professional, Innovational

Wireless Electric Derma Pen BL6

Product description:

- Tightens, lifts and rejuvenates the skin
- Reduces acne and surgical scarring
- Improves wrinkles & fine lines
- Minimises pore size
- Improves stretch marks


Specifications

Type	Derma Rolling System
Function:	Microneedle Therapy System
Handpiece material	Aviation aluminum
Needle material	Titanium & Stainless Steel
Needles	9-12-36-42- Nano
Needle length	0.25- 2.5mm adjustable
Speed control	5 Level
Adapter	US/EU/UK/AU
Input voltage	110-220v
Output voltage	5V/1A
Power mode	Wireless/wired
Rechargeable battery	2pcs
Battery working time	4 hours
Weight	Main body 56g
Warranty	12 months

Standard configuration

1. Host	1 pc
2. Adapter	1 pc
3. Cord	1 pc
4. Cover	1 pc
5. Battery	2 pcs
6. Charger	1 pc
7. Manual	1 pc
8. Tray	1 pc
9. Cartridges	2 pcs 12pins

BEFORE AND AFTER



Wuxi Biomedical Technology Co., Ltd.

www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobile phone: 86-13382885910

is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.

